







Standard Product Range Rods

Diameter  mm	Rod weight Length approx. 1,500 mm g 	Carton contents		Pallet load	
		Number of rods	Weight approx. kg	Number of boxes	Weight approx. kg
3 ± 0.13	24	529	12.5	44	550.0
4 ± 0.13	42	298	12.5	44	550.0
5 ± 0.13	66	183	12.0	44	528.0
6 ± 0.13	95	140	13.2	44	580.8
7 ± 0.13	129	98	12.6	44	554.4
8 ± 0.18	168	80	13.4	44	589.6
9 ± 0.18	213	63	13.4	44	589.6
10 ± 0.18	263	45	11.8	44	519.2
12 ± 0.18	378	35	13.2	44	580.8
14 ± 0.26	515	24	12.4	44	545.6
16 ± 0.26	672	20	13.4	36	482.4
18 ± 0.36	851	20	17.0	27	459.0
20 ± 0.36	1 050	16	16.8	27	453.6
22 ± 0.40	1 271	12	15.3	36	550.8
24 ± 0.40	1 512	12	18.2	27	491.4
26 ± 0.50	1 775	9	16.0	27	432.0
28 ± 0.70	2 059	9	18.5	27	499.5
30 ± 0.70	2 363	6	14.2	36	511.2

Standard length: approx. 1,500 mm

Standard Product Range Capillaries

Outside diameter  mm	Inside diameter  mm	Tube weight Length approx. 1,500 mm g 	Carton contents	
			Number of tubes	Weight approx. kg
4 ± 0.16	0.8 ± 0.08	40	250	10.0
5 ± 0.16	0.4 ± 0.08	65	154	10.0
	0.6 ± 0.08	65	154	10.0
	0.8 ± 0.08	64	156	10.0
6 ± 0.16	1.2 ± 0.08	62	161	10.0
	0.4 ± 0.08	94	104	10.0
	0.8 ± 0.08	93	108	10.0
	1.2 ± 0.08	91	110	10.0
7 ± 0.18	1.7 ± 0.10	87	115	10.0
	2.2 ± 0.10	82	122	10.0
	2.7 ± 0.10	75	133	10.0
	0.8 ± 0.08	127	79	10.0
8 ± 0.18	1.2 ± 0.08	125	80	10.0
	1.7 ± 0.10	121	83	10.0
	2.2 ± 0.10	116	86	10.0
	2.7 ± 0.10	110	91	10.0
	3.0 ± 0.10	105	95	10.0
9 ± 0.18	0.8 ± 0.08	166	60	10.0
	1.2 ± 0.08	164	61	10.0
	1.7 ± 0.10	160	62	10.0
	2.2 ± 0.10	155	64	10.0
	2.7 ± 0.10	149	67	10.0
9 ± 0.18	3.0 ± 0.10	144	69	10.0
	0.8 ± 0.08	211	47	10.0
	1.2 ± 0.08	209	48	10.0
	1.7 ± 0.10	205	49	10.0
9 ± 0.18	2.2 ± 0.10	200	50	10.0
	2.7 ± 0.10	194	52	10.0
	3.0 ± 0.10	189	53	10.0

 Pallet loading capillaries:
number of cartons: 55
weight: approx. 550.0 kg

Standard length: approx. 1,500 mm